



CelCulture® 170L CO₂ Incubators




CelMate® 170L CO₂ Incubators



CelSafe® 170L CO₂ Incubators

Esco CO₂ Incubators

Safe Cradle for Cell Cultivation



Designed in the USA



Table of Contents

About Esco	03
Esco Products Overview	05
Esco CO ₂ Incubators	
• CelCulture®	06
• CelCulture® with UV Lamp	16
• CelCulture® with Copper Interior Chamber	17
• CelCulture® with Integrated Cooling System	18
• CelCulture® with Stainless Steel Exterior	19
• CelCulture® for <i>In Vitro</i> Fertilization	20
• CelCulture® with Water Jacket	23
• CelMate®	25
• CelSafe®	27
Options and Accessories	32
CelCradle™ Bioreactor	36
After Sales Services	39





Welcome to Esco

Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.



Esco represents innovation and forward-thinking designs, which are all coupled with the highest standard quality since 1978. The Esco Group of Companies remains dedicated in delivering innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical and IVF community. With the most extensive product line in the industry, our products have passed a number of international standards and certifications. Esco operates under ISO 9001, ISO 14001 and ISO 13485.

Availability and Accessibility. Headquartered in Singapore, manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, services and marketing subsidiaries are located in 12 major markets including the US, UK, Japan, China and India. Our regional distribution centers are located in Bangladesh, China, Denmark, Germany, Hong Kong, India, Indonesia, Italy, Japan, Lithuania, Malaysia, Philippines, Russia, Singapore, South Africa, South Korea, Russia, Taiwan, Thailand, UAE, United Kingdom, USA, and Vietnam. Because of our worldwide presence, you can be sure that Esco is within your reach.

High Quality, Reliable and Dependable. Our customers are confident that only with the best quality, reliable, and dependable products, can they be sure of the accuracy of their research and procedures. Cross functional teams from Esco Production, R&D, Quality Assurance and Senior Management, are regularly assembled to review and implement areas for improvement.

Esco Cares for Your Safety. Esco focuses on providing safety not just for your samples but also for you and the environment.


Esco Cares for Your Comfort. Comfort of our users is ensured by building ergonomic designs and by reducing the noise levels of the units.

Esco Cares for the Environment. One in every four of Esco's employees is involved in Research and Development and a number of these evaluate new components and/or designs to produce energy efficient equipment. Being GREEN is more than just modifying the parts we use to produce a new energy efficient technology, it also embodies the every aspect of the company.

Customer Service and Support. Our service does not stop once purchase has been made. Esco gives on-time customer service and offers end-user seminars, service training, preventive maintenance, provides educational materials and informative videos.

As Esco takes the opportunity to respond to the world's needs, we aim not just to contribute in the advancement of scientific discoveries but also in making the world a safer, healthier and a better place to live in.

Research and Development



An integral part of our business planning effort is based on managing a robust research and development program in Singapore, China, Europe and the USA, balanced against an investment in service support, training and customer education. Compared to industry averages, Esco invests a significant percentage of annual revenues in research and development. As a result of our investment, and with continuous feedback and idea evaluation among our research, global sales, marketing, purchasing and manufacturing teams; Esco products reflect the best contemporary designs in performance, ergonomics and customer satisfaction.

- Engineers located in technology centers in Singapore, China, Europe and the USA
- Growing patent portfolio
- Core competencies:
 - Embedded system, sensor and software development and integration
 - Containment engineering for biohazards, chemical vapors and hazardous powders
 - Decontamination cycle development
 - Computational fluid dynamics
 - Temperature, humidity, gas and environmental control
 - Imaging systems
 - Wireless and remote monitoring
 - cGMP laboratory design

Production and Quality



Esco's manufacturing advantage stems from our extensive degree of vertical integration, enabled by our world-leading high throughput. All processes, with a few exceptions, are performed in-house. This allows us to achieve quality and reliability that is truly world class.

- Incoming materials inspection and warehousing
- CNC-controlled sheet metal fabrication
- Full range of welding including MIG, TIG, Spot, and Robotic welding
- Environment-friendly powder coating lines
- Electromechanical final product assembly
- Electrical / electronics sub-assembly
- Multi-step electrical and physical performance testing
- Independent quality control at each step in the production cycle
- Microbiology, chemistry, containment test labs
- Instruments calibration laboratory

Esco's focus on quality and timeliness is relentless. Continuous improvement is a mantra. Cross-functional teams from Esco Production, R&D, Quality Assurance, Senior Management, are regularly assembled to review and implement areas for improvement.

Every year, Esco manufacturing site is audited by certification bodies like ISO 9001, ISO 14001, ISO 13485, NSF International, and TUV Nord Germany. Our machines safety, workplace safety, and environmental safety are also audited by different government institutions.

All of the third party audits help our customer to ensure Esco manufacturing sites consistency of producing quality products and continual improvement.

Products and Application

Life Sciences Laboratory Equipment

Sample Preparation

- Class I Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Clean Benches
- Vertical Laminar Flow Clean Benches
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

Sample Cultivation

- CO₂ Incubators, Direct Heat Air-Jacketed
- CO₂ Incubators with Cooling System
- CO₂ Incubators with Stainless Steel Exterior
- CO₂ Incubators (Water-Jacketed)
- Laboratory Shakers

Sample Analysis

PCR Thermal Cyclers

- Conventional Thermal Cyclers
- Real-time PCR Systems

PCR Sample Handling

- Microplate Shakers
- PCR Cabinets

Sample Storage & Sample Protection Solutions

- Ultra-low Temperature Freezers
- Lab Refrigerators and Freezers
- Sample Database Management Software
- Intelligent Remote Monitoring Application Protocol
- Remote Monitoring, Datalogging, Programming Software
- Wireless Monitoring System

Chemical Research

- Ductless Fume Hoods
- Laboratory Fume Hoods
- Fume Hood Airflow Monitors
- Exhaust Blowers
- Powder Weighing Balance Enclosures

General Equipment

Laboratory Thermostatic Products

- Laboratory Oven
- Forced Convection Laboratory Incubator
- Natural Convection Laboratory Incubator
- Refrigerated Incubator

Forensic Sciences

- Evidence Drying Cabinet

Medical / IVF Equipment

Controlled Embryo Handling

- Fertilisafe™ ART Workstation
- AVT-I Anti Vibration Table
- Semi Closed Environment IVF

Safe Embryo Culture

- Miri® Multi room Benchtop Incubator
- CelCulture® CO₂ Incubator
- Mini Miri® Humidified Benchtop Incubator

Innovative Time Lapse Imaging

- Miri® Time-lapse Incubator

Accurate Quality Control

- Miri® GA Gas and Temperature Validation Unit
- Miri® GA Mini Gas Validation Unit

Unique Consumables

- CultureCoin

Healthcare

Esco Pharma Products

Airflow Containment Products

- Pharmacon® Downflow Booths
- Ceiling Laminar Airflow Units
- Laminar Flow Horizontal/Vertical Trolley
- Enterprise™ Laminar Flow Straddle Unites
- Cytoculture™ Cytotoxic Safety Cabinets

Isolation Containment

- Aseptic Containment Isolator (ACTI)
- Weighing and Dispensing Containment Isolator (WDCI)
- General Processing Platform Isolator (GPPI)
- Containment Barrier Isolator (CBI)
- Turbulent Flow Aseptic (Grade A) Isolator (TFAI)
- Isoclean® Healthcare Platform Isolator (HPI)
- Streamline® Compounding Isolators (SCI)
- Technetium Dispensing Isolators
- Blood Cell Labeling Isolators
- Open and Closed Restricted Barrier Access Systems (RABS)

Cross Contamination Facility Integrated Barrier

- BioPass™ Pass Through
- Infinity® Air Shower Pass Box
- Cleanroom Air Shower
- Infinity® Cleanroom Transfer Hatch
- Infinity® Pass Box
- Soft capsule® Soft Wall Cleanroom
- Dynamic Passboxes and Dynamic Floor Label Hatches
- Laminar Flow Storage Cabinet

Ventilation Containment

- Ventilated Balance Enclosure

VacciXcell Products

Bioreactors and Fermenters

- CelCradle™
- TideCell®
- VacciXcell™ Hybrid Bioreactor

Cell Culture Monitoring, Media and Consumables

- Super Plus™
- Plus™ Vero
- Plus™ MDCK
- Plus™ MDCK II
- BioNOC™ II macrocarriers
- GlucCell® Glucose Monitoring System
- CVD Kit

Filling Line Equipment

- Filling Line Isolators
- cRabs (close restricted access barriers)
- oRabs (open restricted access barriers)

Integrated Solutions

- Cell Processing Isolator
- Cell Processing Center

TaPestle Rx Products and Services

PRODUCTS

Pharmacy Automation and Compounding Supply

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (Class II BSC, VBE, LFC)
- Radiopharmacy Hoods and Isolators
- Aseptic Filling Systems

Healthcare and Laboratory Construction Components

- Prefabricated Walls (Airecell®)
- Prefabricated Containerized Facility (Prefab™)
- Series Ceiling Systems
- Hygienic/Hermetic Door Systems
- Surgical Scrub Sinks
- Vinyl Tiles and Epoxy
- Laboratory Fit-outs
 - Worktops
 - Frames
 - Specialty Storage cabinets
 - Service Spines & Reagent Shelving

SERVICES

- Conceptualization
- Planning
- Procurement
- Installation

FACILITY DESIGNS

- Process Architecture
- Biocontainment/Biosafety
- Pharmacy Compounding/Nuclear Medicine
- Cleanroom, Vaccine and Cell Processing
- Laboratory
- Containerized Facility
- ART/IVF
- Cold Chain



CelCulture® CO₂ Incubators

INTRODUCTION

CO₂ incubators are widely used in scientific research to grow and maintain cell cultures. Typical fields of application include tissue engineering, *in vitro* fertilization, neuroscience, cancer research and other mammalian cell research applications.

Sleek, reliable and intuitive, Esco CelCulture® CO₂ incubators provide complete sample protection that brings your scientific dreams one step closer to reality.

KEY FEATURES


CelCulture® CO₂ INCUBATORS *Cradle for Beautiful Cells*



CelCulture® CO₂ Incubators
available in 3 sizes, 50 L, 170 L, and 240 L.

ULPA FILTER*

- 99.999% efficient, superior to conventional HEPA filters
- Filters air continuously
- Chamber returns to ISO Class 5 cleanliness in 11 minutes upon door closing to prevent contamination



SHELVING

- Perforated shelving to improve uniformity
- Anti-tip
- Stainless steel
- Built-in grip
- Dismantles without tools for easy cleaning

DIRECT HEAT & AIR JACKET

- Fast and uniform heating
- Rapid temperature recovery
- Air jacket improves chamber stability



DUCT WORK

- Directs air flow for rapid recovery and excellent uniformity
- Easily removed for cleaning



WATER PAN

- Precisely heated by base heater to provide high humidity
- Gentle airflow over water surface accelerates humidity recovery



ROUNDED CORNERS

- Seamless design
- Facilitates easier cleaning

O₂ SENSOR

for suppressed O₂ model

- Long life
- Stable output signal
- Not affected by CO₂ concentration inside the chamber

DOOR SWITCH

Automatically turns off the blower and gas supply when the door is opened

TOP COVER

Provides quick access to electrical panel components

CO₂ SENSOR

- IR Sensor
- Single-beam, dual-wavelength IR sensor is drift-free
- Auto-zeroing

SMARTSENSE™ MICROCONTROLLER INTERFACE

Intuitive, fully-equipped control and monitoring system

BLOWER*

Gentle airflow in chamber improves recovery and uniformity

OUTER DOOR

- Reversible
- Heated to prevent condensation

SAMPLE PORT

Allows direct measurement of chamber atmosphere such as temperature and CO₂ concentration

INNER GLASS DOOR

For observing sample cells inside the chamber during operation

DOOR LATCH

To lock / unlock the glass door

LEVELING FEET

Easily adjustable


QUALITY ESCO CONSTRUCTION

- Electrogalvanized steel with white oven-baked epoxy-polyester antimicrobial powder-coated finish.
- External surfaces are powder coated with Esco **ISOCIDE™** to eliminate 99.9% of surface bacteria within 24 hours of exposure.
- Ensures a healthier, safer and cleaner lab environment.


**Not available in 50 L Model (including the top and bottom plenum)*

VIVOCELL™ PRECISE PARAMETER CONTROL

BEST UNIFORMITY AND CONTROL AMONG THE COMPETITION




Different lines represent different sensor positions inside the chamber. Esco CelCulture® has uniformity variance of less than $\pm 0.5^\circ\text{C}$ which means all the samples are evenly heated.*



Minimal fluctuation ($\pm 0.2^\circ\text{C}$) ensures temperature stability.*

FAST CO₂, TEMPERATURE AND HUMIDITY RECOVERY WITHOUT OVERSHOOT





Precisely-tuned sensor and software result in fast recovery of CO₂ without overshoot. This ensures uniform CO₂ levels even with frequent incubator door opening.

Recovery of both temperature and humidity is twice as fast as conventional incubators.

- Company A's model: overshoot.
- Company B's model: slow recovery.
- Esco CelCulture®: fast recovery, no overshoot.



DIRECT HEAT AND AIR JACKET



- Direct heating enables rapid temperature recovery while air jacket provides isolation against ambient temperature fluctuations.
- Precise heating in the chamber is achieved by using 8 heaters located in 3 zones. The 3 zones are intelligently controlled by the microcontroller for best temperature uniformity and minimal fluctuation.

- The main heater provides precise temperature control.
- The bottom heater warms the water pan and provides humidity.
- The outer door heater prevents condensation on glass door and facilitates temperature recovery.

VENTIFLOW™ FORCED CONVECTION




- No disturbance to cell culture.
- Blower automatically stops when door is opened to minimize mixing of chamber and room air.
- Accelerates recovery of chamber air to ISO Class 5 Cleanliness after door closing to prevent contamination.
- Improves CO₂, humidity and temperature uniformity.
- Filtered air circulates across water pan to accelerate humidifying process.

* Units were factory-tested under controlled environmental conditions per Esco method. Esco does not guarantee identical results in the field under differing conditions. Original report available upon request. Model used in the test is CCL-170B-8.

ROBUST CONTAMINATION CONTROL


STERISAFE™ ULPA FILTRATION SYSTEM



- Chamber air is continuously filtered by ULPA filters to keep the chamber at ISO Class 5 cleanliness. This ensures that all contaminants from both room air and chamber air are filtered, thus only clean air is recirculated.
- ULPA filters operate at 99.999% efficiency, superior to conventional HEPA filters which are 99.99% efficient.
- Chamber achieves ISO Class 5 cleanliness 11 minutes after door closing.*

* Units were factory-tested under controlled environmental conditions per Esco method. Esco does not guarantee identical results in the field under differing conditions. Original report available upon request. Model used in the test was CCL-170B-8.

VALIDATED SWIFTCON™ OVERNIGHT DECONTAMINATION CYCLE



Microorganisms	Before Decon	After Decon
<i>Bacillus atrophaeus</i>	1.59 x 10 ⁶	0
<i>Aspergillus brasiliensis</i>	1.52 x 10 ⁴	0
<i>Pseudomonas aeruginosa</i>	2.38 x 10 ⁶	0
<i>Staphylococcus epidermis</i>	2.33 x 10 ⁶	0
<i>Escherichia coli</i>	1.57 x 10 ⁶	0
<i>Staphylococcus aureus</i>	5.72 x 10 ⁶	0
<i>Enterobacter faecalis</i>	2.15 x 10 ⁶	0

- The automated SwiftCon™ 90°C moist heat decontamination cycle has been proven effective in deactivating normally resistant fungi, bacterial spores and vegetative cells by the Health Protection Agency (HPA) in UK.
- Full decontamination cycle completes within 15 hours.
- Chamber is cool and dry at the end of the cycle. No further wipe down is needed.*

- Independently proven to be as effective as high temperature decontamination.
- Lower temperature causes less damage to electronic components, therefore prolongs the life span of the incubator.

*Not applicable to CCL-50L unit since it has no decon pump to dry the chamber and condensation will normally occur in the chamber after the cycle. Further wipe down is therefore required after the cycle is done.


GAS INJECTION LINES ARE FILTERED



- All gas injection lines are filtered via 0.2 micron inlet filters to remove impurities and contaminants before being injected into the chamber.
- Inlet filters are field-replaceable and are located external to the incubator.

CONTROLLER TYPE

USER - FRIENDLY SOFTWARE INTERFACE



1. **Start / stop decontamination cycle**
2. **HEAT LED**
Lights up when heat is applied to chamber
3. **°C is lit when displaying the temperature**
4. **% RH is lit when displaying the humidity level**
5. **Enter menu / go back to previous menu**
6. **Scroll up / increase value**
7. **ALARMS LED**
Will blink when errors and warnings occur
8. **Mute alarms**
9. **INJECT LED**
Lights up when gas is injected
10. **% O₂ is lit when displaying the O₂ concentration**
11. **% CO₂ is lit when displaying the CO₂ concentration**
12. **Confirm value / enter a menu**
13. **Scroll down / decrease value**

• Comprehensive, user-configurable alarms:

- Temperature
- CO₂
- Humidity (if installed)
- O₂ (if installed)

• CelAlert™ alarm system reminds user to replace parts.

```
NEW CO2 TANK
CO2 TANK NET WEIGHT
IN KG: 0031
NEW CO2 TANK IS SET
```

In addition to CO₂ tank low alarm, CelAlert™ has CO₂ tank depletion reminder that automatically calculates how much CO₂ gas is left in the tank and alerts user to replace the tank one week before the gas is depleted. This gives the user some buffer time to place orders for new tanks.

```
SET ULPA REMINDER
ULPA REMINDER
TIME <MONTHS> :12
REMINDER IS SET
```

ULPA reminder will alert user to replace ULPA filter.

• Intelligent data and event logger records all incubator parameters for on-screen recall. A 2 MB built-in flash memory guarantees long-term storage of data.

```
MAIN MENU
CALIBRATION
ADMIN
>DATALOGGING
```

```
DATALOGGING
>SHOW TEMP DATA LOG
SHOW %CO2 DATA LOG
DATA LOG PERIOD
```

```
012016 0724 36.8°C
012016 0719 37.0°C
012016 0714 37.1°C
012016 0709 37.3°C
```

• Diagnostic interface and online quick help provide comprehensive solutions to frequently encountered problems.

Voyager®

Remote Monitoring, Datalogging, Programming Software

Esco Voyager® is a PC-based software package developed for the remote monitoring, datalogging, and programming / device configuration of Esco thermostatic products.

A centralized monitoring and control system for the laboratory, Esco Voyager® provides extra protection for you and your samples.



Voyager® interfaces with individual Esco equipment over RS485 using the EscoBUS communications protocol. Multiple equipment may be interfaced to a single PC.

Compatible Equipment

- Lexicon® - Ultra-low Temperature Freezer
- CelCulture® - CO₂ Incubator (CCL)
- CelMate® - CO₂ Incubator (CLM)
- CelSafe® - CO₂ Incubator (CLS)
- Isotherm® - Forced Convection Oven (OFA)
- Isotherm® - Forced Convection Incubator (IFA)
- Isotherm® - Refrigerated Incubator (IFC)
- Isotherm® - Natural Convection Incubator (INA)



REAR PANEL



1 Cooling Fan*

Prevents the electrical panel from overheating.

**Available in 170 L and 240 L models only*



2 Power Supply Inlet

Connects the incubator unit to the power source.



6 N₂ Gas Supply Inlet**

Connects the N₂ gas supply to the incubator. Inlet pressure requirement is 15 psi.

***Applicable only to Suppressed O₂ models*



7 CO₂ Gas Supply Inlet

Connects the CO₂ gas supply to the incubator. Inlet pressure requirement is 15 psi.




3 RS485 Communication Port

Provides serial communication port for PC. It can be daisy-chained from one product to another and can also be connected to a PC



8 Gas Inlet Filter

Provided to remove any contaminants from the gas supply.



4 Analog Port (Optional)


Allows the incubator to output analog signals representing temperature, CO₂/O₂** concentration and relative humidity, depending on the options available in the incubator. This allows the incubator to be connected to an in-house data acquisition or alarm system.



5 Alarm Contact

A set of relay contacts located on the rear panel of the unit is provided to monitor temperature, humidity, O₂** or CO₂ alarms. These can be connected to a remote alarm system.

***Applicable only to Suppressed O₂ models*



9 Access Port

Allows cables, hoses or additional sensors to be routed into the work space. A rubber stopper with controlled leak is installed as standard configuration and is part of standard accessories.

CelCulture® CO₂ INCUBATOR SENSORS



IR SENSOR

An Infrared (IR) sensor is a versatile instrument for measuring CO₂ level inside the incubator. The CARBOCAP® sensor is silicon-based and its operation is based on the NDIR Single-Beam Dual-Wavelength principle.

IR-based sensors are not affected by water vapor, dust or most chemicals. The single-beam dual-wavelength technology (one reference and one measurement) ensures a drift-free sensor that does not require calibration by the user.

Operating principle

The light source is positioned to shine at the IR detector so that the light travels a fixed distance to the detector, where the intensity of the light is measured. A Fabry-Perot Interferometer (FPI) is positioned just in front of the IR detector. The FPI is a tunable filter which allows only certain wavelengths of light to pass through to the detector.

Carbon dioxide absorbs certain wavelengths of light and not others, so the FPI is designed to pass light at a CO₂ absorption wavelength (4.26 µm) and a nearby, non-absorbing wavelength.

When the sensor is operating, the FPI is regularly tuned back and forth between the two wavelengths. At the CO₂ absorption wavelength, the intensity of detected light is reduced in proportion to the concentration of CO₂ in the optical path. The light intensity measured at the non-absorbing wavelength serves as a baseline for comparison.

Operating Conditions:

%CO₂ detection range: 0 to 20% CO₂ Concentration

%RH operating range: Not affected by Humidity

Temperature range: -20°C to +60°C



O₂ SENSOR

Figaro's O₂ sensor is a unique galvanic cell type oxygen sensor. Its most notable features include long life expectancy and excellent chemical durability, and it is not influenced by CO₂. This O₂ sensor ideally meets the ever-increasing demand for oxygen monitoring in various fields including Biochemistry, Limnology, Medicine, Soil Respiration and in combustion gas monitoring.


Operating Conditions:

%O₂ detection range: 1 to 20.7% O₂ Concentration

%RH operating range: 10% to 90% Relative Humidity

Temperature range: 5°C to 40°C

Technical drawing showing the location of the IR sensor and O₂ sensor in relation to the other chamber components of the CelCulture® CO₂ Incubator.



ORDERING INFORMATION

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-050B-8	2170034	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz, (Without Decon Pump)
CCL-170B-8	2170002	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz,
CCL-170B-8-NF	2170068	CelCulture® Incubator 170 L IR Sensor, CO ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz, (No ULPA Filter)
CCL-240B-8	2170058	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-240B-8-NF	2170069	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz, (No ULPA Filter)
CCL-050B-9	2170054	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9	2170004	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-170B-9-NF	2170075	CelCulture® Incubator 170 L IR Sensor, CO ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz, (No ULPA Filter)
CCL-240B-9	2170060	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-240B-9-NF	2170079	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz, (No ULPA Filter)

SURPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER


MODELS	ITEM CODE	DESCRIPTION
CCL-050T-8	2170055	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8	2170047	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-170T-8-NF	2170070	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz (No ULPA Filter)
CCL-240T-8	2170061	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-240T-8-NF	2170071	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz (No ULPA Filter)
CCL-050T-9	2170056	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9	2170048	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-170T-9-NF	2170076	CelCulture® Incubator 170L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ (No ULPA Filter)
CCL-240T-9	2170062	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-240T-9-NF	2170080	CelCulture® Incubator 240L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ (No ULPA Filter)

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER WITH FLAT DOOR DESIGN

MODELS	ITEM CODE	DESCRIPTION
CCL-050B-8-FD	2170150	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-8-FD	2170117	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz
CCL-170B-8-NF-FD	2170243	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz, (No ULPA Filter)
CCL-240B-8-FD	2170123	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz
CCL-240B-8-NF-FD	2170244	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz, (No ULPA Filter)
CCL-050B-9-FD	2170191	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9-FD	2170120	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz
CCL-240B-9-FD	2170126	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER WITH FLAT DOOR DESIGN

MODELS	ITEM CODE	DESCRIPTION
CCL-050T-8-FD	2170149	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8-FD	2170118	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz
CCL-170T-8-NF-FD	2170247	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz (No ULPA Filter)
CCL-240T-8-FD	2170125	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz
CCL-240T-8-NF-FD	2170249	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Flat Door, 230 VAC, 50/60 Hz (No ULPA Filter)
CCL-050T-9-FD	2170245	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9-FD	2170148	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz
CCL-240T-9-FD	2170127	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, Flat Door, 115 VAC, 50/60 Hz



CelCulture® CO₂ Incubator with UV Lamp*

INTRODUCTION

Ozone-free UV lamp decontaminates the water in the humidity pan and the circulating air that flows across the water surface. The UV lamp is isolated from the inner chamber by the plenum cover, ensuring the safety of cell cultures from the UV decontamination.

WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-170B-8-UV	2170031	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170B-9-UV	2170121	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-170T-8-UV	2170043	CelCulture® Incubator 170 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170T-9-UV	2170044	CelCulture® Incubator 170 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240B-8-UV	2170094	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240B-9-UV	2170096	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240T-8-UV	2170097	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240T-9-UV	2170098	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz

WITH COPPER CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-170B-8-Cu-UV	2170099	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170B-9-Cu-UV	2170285	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240B-8-Cu-UV	2170154	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240B-9-Cu-UV	2170286	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-170T-8-Cu-UV	2170131	CelCulture® Incubator 170 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170T-9-Cu-UV	2170287	CelCulture® Incubator 170 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240T-8-Cu-UV	2170155	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240T-9-Cu-UV	2170288	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz

WITH PELTIER COOLING MODE

MODELS	ITEM CODE	DESCRIPTION
CCL-170B-8-P-UV	2170204	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170B-9-P-UV	2170205	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-170T-8-P-UV	2170206	CelCulture® Incubator 170 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-170T-9-P-UV	2170207	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240B-7-P-UV	2170289	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 100 VAC, 50/60 Hz
CCL-240B-8-P-UV	2170208	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240B-9-P-UV	2170209	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz
CCL-240T-7-P-UV	2170262	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 100 VAC, 50/60 Hz
CCL-240T-8-P-UV	2170210	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 220 - 240 VAC, 50/60 Hz
CCL-240T-9-P-UV	2170211	CelCulture® Incubator 240 L, IR Sensor, CO ₂ , O ₂ Control, 90°C Moist Heat Decon, 110 - 130 VAC, 50/60 Hz

*Not for sale in North America, Europe, Korea, Japan, Australia and New Zealand



CelCulture®

CO₂ Incubator with Copper Interior Chamber

Pure solid copper interior offers additional protection for your precious samples.

MAXIMUM CONTAMINATION CONTROL

Copper has been known for millennia to have antimicrobial properties. Copper can inhibit the growth of common culture microbial contaminants such as:

- *Escherichia coli*
- *Staphylococcus aureus*
- viruses

Other contamination control methods include:

- ✓ ULPA filter with 99.999% efficiency*
- ✓ 90°C Moist Heat Decontamination Cycle (HPA-Validated)
- ✓ 0.2 micron inlet filter for gas inputs
- ✓ ISOCIDE™ antimicrobial powder coating

**Not available in 50 L model*

IR SENSOR MODEL WITH 100% COPPER CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-050B-8-Cu	2170081	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-8-Cu	2170083	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-240B-8-Cu	2170085	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-050B-9-Cu	2170082	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9-Cu	2170084	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-240B-9-Cu	2170086	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH COPPER CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-050T-8-Cu	2170087	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8-Cu	2170089	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-240T-8-Cu	2170111	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz
CCL-050T-9-Cu	2170088	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9-Cu	2170090	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz
CCL-240T-9-Cu	2170092	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz



CelCulture® CO₂ Incubator with Integrated Cooling System

INTRODUCTION

Esco CelCulture® CO₂ Incubator with Integrated Cooling System provides solution for highly specialized applications.

The integrated cooling system allows studies of samples that requires temperature at/or below ambient temperature.

KEY FEATURES

WIDER TEMPERATURE RANGE

Temperature range of 12°C below ambient to 60°C above ambient means wider range of applications.

HIGHLY EFFICIENT, ENVIRONMENT FRIENDLY Peltier COOLING SYSTEM

This provides precise heating and cooling inside the chamber making sure that your samples are safe from temperature changes.

COMPLETE CONTAMINATION CONTROL METHODS

- 90 °C validated moist heat decontamination cycle
- ULPA filter
- ISOCIDE™ anti-microbial coating
- 0.2 micron inlet filter for gas inputs

IR SENSOR MODEL WITH INTEGRATED COOLING SYSTEM

MODELS	ITEM CODE	DESCRIPTION
CCL-170B-8-P	2170101	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Peltier System, 230 VAC, 50/60 Hz
CCL-170B-9-P	2170115	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Peltier System, 115 VAC, 50/60 Hz
CCL-240B-8-P	2170116	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Peltier System, 230 VAC, 50/60 Hz
CCL-240B-9-P	2170266	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Peltier System, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH INTEGRATED COOLING SYSTEM

MODELS	ITEM CODE	DESCRIPTION
CCL-170T-8-P	2170112	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Peltier System, 230 VAC, 50/60 Hz
CCL-170T-9-P	2170153	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Peltier System, 115 VAC, 50/60 Hz
CCL-240T-8-P	2170267	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Peltier System, 230 VAC, 50/60 Hz
CCL-240T-9-P	2170268	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Peltier System, 115 VAC, 50/60 Hz



CelCulture®

CO₂ Incubator with Stainless Steel Exterior Cabinet

- Corrosion-resistant Surface
- Meets Pharmaceutical and Clinical Laboratory Requirements

IR SENSOR MODEL WITH STAINLESS STEEL EXTERIOR CABINET

MODELS	ITEM CODE	DESCRIPTION
CCL-050B-8-SS	2170128	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-8-SS	2170065	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240B-8-SS	2170137	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-050B-9-SS	2170176	CelCulture® Incubator 50 L, IR sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9-SS	2170177	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz
CCL-240B-9-SS	2170140	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL EXTERIOR CABINET

MODELS	ITEM CODE	DESCRIPTION
CCL-050T-8-SS	2170171	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8-SS	2170129	CelCulture® Incubator 170 L IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240T-8-SS	2170138	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-050T-9-SS	2170178	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9-SS	2170179	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz
CCL-240T-9-SS	2170141	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz



CelCulture®

CO₂ Incubator for *In Vitro* Fertilization

The CO₂ Incubator has a vital role in providing an optimal environment in embryo development during IVF and other ART procedures.

Sleek, reliable and intuitive, Esco CelCulture® CO₂ incubators for IVF are packed with outstanding features such as rapid parameter recovery, ISO Class 5 Cleanliness, ISOCIDE™ antimicrobial coating, standard Inner Door Kit feature that reduces contamination risk, and other accessories for specialized applications.



CCL-050_--IVF

KEY FEATURES

- Wider temperature range, from ambient +3°C temperature to 60°C.
- Complete contamination control methods to protect your precious samples.
- All gas inputs are filtered via 0.2 micron inlet filter and has ULPA filtration system.
- 90°C moist heat decontamination cycle, validated by HPA-UK.

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODEL	ITEM CODE	DESCRIPTION
CCL-050B-8-IVF	2170257	Celculture® Incubator, 50L, IR sensor, CO ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 50L (2 Glass Doors with Latches), Factory Installed, 230VAC, 50/60 Hz
CCL-170B-8-IVF	2170272	CelCulture® Incubator 170L IR Sensor, CO ₂ Control ULPA, Moist Heat Decon, with Sealed Inner Door Kit for 170L (4 Glass Doors with Latches), Factory Installed, 230VAC 50/60Hz
CCL-240B-8-IVF	2170278	CelCulture® Incubator 240L IR Sensor CO ₂ Control, ULPA, Moist Heat Decon, with Sealed Inner Door Kit for 240L (4/6 Glass Doors with Latches), Factory Installed, 230VAC 50/60Hz
CCL-050B-9-IVF	2170258	Celculture® Incubator, 50L, IR sensor, CO ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 50L (2 Glass Doors with Latches), Factory Installed, 115VAC, 50/60 Hz
CCL-170B-9-IVF	2170273	CelCulture® Incubator 170L IR Sensor, CO ₂ Control ULPA, Moist Heat Decon, with Sealed Inner Door Kit for 170L (4 Glass Doors with Latches), Factory Installed, 115VAC 50/60Hz
CCL-240B-9-IVF	2170279	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon with Sealed Inner Door Kit for 240L (4/6 Glass Doors with Latches), Factory Installed, 115VAC 50/60Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER

MODEL	ITEM CODE	DESCRIPTION
CCL-050T-8-IVF	2170260	Celculture® Incubator, 50L, IR sensor, CO ₂ & O ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 50L (2 Glass Doors with Latches), Factory Installed, 230VAC, 50/60 Hz
CCL-170T-8-IVF	2170275	CelCulture® Incubator 170L IR Sensor, CO ₂ & O ₂ Control ULPA, Moist Heat Decon, with Sealed Inner Door Kit for 170L (4 Glass Doors with Latches), Factory Installed, 230VAC 50/60Hz
CCL-240T-8-IVF	2170281	Celculture® Incubator, 240L, IR sensor, CO ₂ & O ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 240L (4/6 Glass Doors with Latches), Factory Installed, 230VAC 50/60Hz
CCL-050T-9-IVF	2170261	Celculture® Incubator, 50L, IR sensor, CO ₂ & O ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 50L (2 Glass Doors with Latches), Factory Installed, 115VAC, 50/60 Hz
CCL-170T-9-IVF	2170276	CelCulture® Incubator 170L IR Sensor, CO ₂ & O ₂ Control ULPA, Moist Heat Decon, with Sealed Inner Door Kit for 170L (4 Glass Doors with Latches), Factory Installed, 115VAC 50/60Hz
CCL-240T-9-IVF	2170282	Celculture® Incubator, 240L, IR sensor, CO ₂ & O ₂ control, Moist Heat Decon, with Sealed Inner Door Kit for 240L (4/6 Glass Doors with Latches), Factory Installed, 115VAC 50/60Hz




TECHNICAL SPECIFICATIONS

Front view




Side view

Rear view




MODEL 50 L



MODEL 170 L



MODEL 240 L



1. Control panel
2. On / off switch
3. Blower fan
4. ULPA filter
5. Sensors

6. Access port
7. Adjustable shelves
8. Humidity pan
9. N₂ gas supply
10. CO₂ gas supply

11. Alarm contact
12. Analog output
13. RS485
14. Power supply inlet
15. Cooling fan

GENERAL SPECIFICATIONS

CEL CULTURE® CO₂ INCUBATORS

CCL-050_-

CCL-170_-

CCL-240_-

TEMPERATURE

Temperature Control Method	Direct heat & Air Jacket using Microcontroller PI		
Ambient Temperature Range	18 to 34°C (64 to 93 °F)		
Temperature Control Range, °C	Ambient +3 to 60		
Temperature Uniformity*, °C	< ±0.5		
Temperature Accuracy, °C	< ±0.1		
Recovery Time** (after 1 min. door opening, 97% from initial value)	≤6 minutes	≤5 minutes	≤5 minutes

CO₂

CO ₂ Control System	Microcontroller PI		
CO ₂ Control Range, % CO ₂	0-20		
CO ₂ Accuracy, % CO ₂	±0.1		
CO ₂ Sensor	Infrared (IR) Sensor		
CO ₂ Recovery Time*** (after 1 min. door opening)	Standard Unit: ≤5 minutes Suppressed O ₂ model: ≤8 minutes	Standard Unit: ≤4 minutes Suppressed O ₂ model: ≤5 minutes	Standard Unit: ≤4 minutes Suppressed O ₂ model: ≤5 minutes

O₂ (FOR SUPPRESSED O₂ MODEL)

O ₂ Control System	Microcontroller PI		
O ₂ Control Range, % O ₂	1-20.7		
O ₂ Accuracy, % O ₂	± 0.1		
O ₂ Sensor	Galvanic Cell Type		
O ₂ Recovery Time (after 1 min. door opening)	At 5.0% O ₂ by volume: 8 minutes	At 5.0% O ₂ by volume: 10 minutes	At 5.0% O ₂ by volume: 12 minutes
	At 1.0% O ₂ by volume: 24 minutes	At 1.0% O ₂ by volume: 24 minutes	At 1.0% O ₂ by volume: 24 minutes

HUMIDITY

Humidification Method	Humidity pan		
Humidity Range****, % RH	Standard Unit: Up to 95% Suppressed O ₂ model: Up to 91%		

PHYSICAL CONSTRUCTION

Interior Volume	50 L (1.8 ft ³)		
External Dimensions (W x D x H)	500 x 500 x 655 mm (19.7" x 19.7" x 25.8")	660 x 660 x 900 mm (26.0" x 26.0" x 35.4")	750 x 770 x 900 mm (29.5" x 30.3" x 35.4")
Internal Dimensions (W x D x H)	345 x 375 x 388 mm (13.6" x 14.8" x 15.3")	505 x 535 x 633 mm (19.9" x 21.1" x 24.9")	595 x 640 x 633 mm (23.4" x 25.2" x 24.9")
Chamber Construction	Main Body	Electrogalvanized steel with ISOCIDE™ antimicrobial coating	
	Interior Material	Stainless steel, type 304	
	Number of Shelves	3	4
	Maximum Number of Shelves	4	7
	Shelves Area	300 x 335 mm (11.8" x 13.2")	465 x 470 mm (18.3" x 18.5")
	Maximum Load per Shelf	4 kg/shelf (8.8 lbs/shelf)	11 kg/shelf (24.3 lbs/shelf)
Electrical Configuration 110-130 VAC, 50/60 Hz	Nominal Power at 37°C	10.6 W	42.9 W
	Maximum Power Consumption	675 W	1184.3 W
	Full Load Amps	5.3 A	9.2 A
Electrical Configuration 220-240 VAC, 50/60 Hz	Nominal Power at 37°C	12.5 W	46.2 W
	Maximum Power Consumption	598.8 W	1008.9 W
	Full Load Amps	2.5 A	4.2 A
Shipping Weight	70 kg (154.3 lbs)	120 kg (264.6 lbs)	155 kg (341.7 lbs)
Shipping Dimensions (W x D x H)	660 x 650 x 900 mm (26.0" x 25.6" x 35.4")	850 x 720 x 1120 mm (33.5" x 28.3" x 44.1")	850 x 850 x 1120 mm (33.5" x 33.5" x 44.1")
Shipping Volume	0.39 m ³ (13.7 ft ³)	0.70 m ³ (24.85 ft ³)	0.79 m ³ (28.03 ft ³)

CONTAMINATION CONTROL

Contamination Control Methods	1) Main body is electrogalvanized steel with ISOCIDE™ antimicrobial coating; 2) 90°C Moist Heat Decontamination Cycle (HPA Validated); 3) 0.2 micron inlet filter for gas inputs; 4) ULPA filter*****		
-------------------------------	--	--	--

*Data recorded under optimum factory setting conditions.

**For temperature not exceeding 37°C .

***For CO₂ not exceeding 5.2%.

**** Esco does not guarantee condensation-free chamber at humidity level higher than 90%.

*****Not available for 50 L.



CelCulture® Water-Jacketed CO₂ Incubators

INTRODUCTION

Esco CelCulture® Water-Jacketed CO₂ Incubator provides a very stable environment to grow and maintain cell cultures.

Designed to protect your cell cultures from sudden power outage and extreme ambient temperature fluctuations, water-jacketed CO₂ incubators maintain constant temperature through the precisely heated water circulating in the jacket walls that surround the chamber.

KEY FEATURES

MORE STABLE TEMPERATURE CONTROL

- Faster temperature recovery times after power outage and door openings.
- Better temperature uniformity

INCREASED SECURITY

Hold a set temperature inside the chamber much longer than air-jacketed units in the event of power failure.

COMPLETE CONTAMINATION CONTROL METHODS

- ULPA filter (for 170/240 L models only)
- ISOCIDE™ antimicrobial coating.
- 0.2 micron inlet filter
- 90°C Moist Heat Decon Cycle (water in the external chamber needs to be drained first)

CONVENIENCE

The unit is equipped with an inlet port to supply water and a drain valve to facilitate faster draining of water when manually cleaning, prior to a decontamination cycle or before transporting the equipment.

EASY MONITORING

Water level can be check via the water level sensor.

WATER-JACKETED, IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-050B-8-WJ	2170156	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-8-WJ	2170103	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz
CCL-240B-8-WJ	2170135	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz
CCL-050B-9-WJ	2170162	CelCulture® Incubator 50 L, IR sensor, CO ₂ Control, Moist Heat Decon, Water-Jacketed, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9-WJ	2170110	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 115 VAC, 50/60 Hz
CCL-240B-9-WJ	2170146	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, Water Jacketed, 115 VAC, 50/60 Hz

WATER-JACKETED, SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CCL-050T-8-WJ	2170132	CelCulture® Incubator 50 L, IR Sensor, CO ₂ and O ₂ Control, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8-WJ	2170157	CelCulture® Incubator 170 L, IR Sensor, CO ₂ and O ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz
CCL-240T-8-WJ	2170144	CelCulture® Incubator 240 L, IR Sensor, CO ₂ and O ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 230 VAC, 50/60 Hz
CCL-050T-9-WJ	2170163	CelCulture® Incubator 50 L, IR sensor, CO ₂ and O ₂ Control, Moist Heat Decon, Water-Jacketed, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9-WJ	2170164	CelCulture® Incubator 170 L, IR Sensor, CO ₂ and O ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 115 VAC, 50/60 Hz
CCL-240T-9-WJ	2170147	CelCulture® Incubator 240 L, IR Sensor, CO ₂ and O ₂ Control, ULPA, Moist Heat Decon, Water-Jacketed, 115 VAC, 50/60 Hz

GENERAL SPECIFICATION

Models		CCL-50-_-_WJ	CCL-170-_-_WJ	CCL-240-_-_WJ
Temperature				
Ambient Temperature Range		18 to 34°C (64 to 93°F)		
Temperature Control Method		Direct Heat and Water Jacketed using Microcontroller PI		
Temperature Range, °C		Ambient +3 to 60		
Temperature Uniformity,* °C		<± 0.5		
Temperature Recovery time** (after 1 min. door opening, 97% from initial value)		≤4 minutes	≤5 minutes	≤5 minutes
Start up time (25°C amb. to 37.0°C)		60 minutes	60 minutes	80 minutes
Power Off Temperature Drop Rate: (1 hour / 10 hours)		1.8°C / 7.0°C	0.6°C / 6.2°C	1.2°C / 5.5°C
CO ₂				
CO ₂ Control System		Microcontroller PI		
CO ₂ Range, %CO ₂		0-20		
CO ₂ Accuracy, %CO ₂		± 0.1		
CO ₂ Sensor		Infrared (IR) Sensor		
CO ₂ Recovery Time*** (after 1 min. door opening)		Standard Unit: ≤6 minutes Suppressed O ₂ model: ≤8 minutes	Standard Unit: ≤4 minutes Suppressed O ₂ model: ≤5 minutes	Standard Unit: ≤5 minutes Suppressed O ₂ model: ≤5 minutes
O ₂ for Supressed O ₂ Model				
O ₂ Control System		Microcontroller PI		
O ₂ Range, %O ₂		1 - 20.7		
O ₂ Accuracy, %O ₂		± 0.1		
O ₂ Sensor		Galvanic Cell Type		
O ₂ Recovery Time*** (after 1 min. door opening)		At 5.0% O ₂ volume: 6 minutes At 1.0% O ₂ volume: 10 minutes	At 5.0% O ₂ volume: 10 minutes At 1.0% O ₂ volume: 20 minutes	At 5.0% O ₂ volume: 12 minutes At 1.0% O ₂ volume: 24 minutes
Humidity				
Humidification Method		Humidity pan		
Humidity Range****, % RH		Standard Unit: Up to 95% Suppressed O ₂ model: Up to 91%		
Physical Parameters				
Interior Volume		50 L (1.8 ft³)	170 L (6.0 ft³)	240 L (8.5 ft³)
Internal Dimensions (W x D x H)		345 x 375 x 388 mm (13.6" x 14.8" x 15.3")	505 x 535 x 633 mm (19.9" x 21.1" x 24.9")	595 x 640 x 633 mm (23.4" x 25.2" x 24.9")
External Dimension (W x D x H)		500 x 525 x 740 mm (19.7" x 20.7" x 29.1")	660 x 680 x 980 mm (26.0" x 26.8" x 38.6")	750 x 795 x 980 mm (29.5" x 31.3" x 38.6")
Water Jacket Volume		12 L	28 L	42 L
Net Weight		75 kg (165 lbs) (no water)	102 kg (225 lbs) (no water)	170 kg (374 lbs) (no water)
Shipping Weight		90 kg (198 lbs)	118 kg (260 lbs)	185 kg (407 lbs)
Shipping Dimension (W x D x H)		660 x 650 x 900 mm (26.0" x 25.6" x 35.4")	850 x 770 x 1200 mm (33.5" x 30.3" x 47.2")	850 x 895 x 1200 mm (33.5" x 35.2" x 47.2")
Number of Shelves		3	4	
Maximum Number of Shelves		4	7	
Shelves Area (W x D)		300 x 335 mm (11.8" x 13.2")	520 x 475 mm (20.5" x 18.7")	550 x 560 mm (21.7" x 22.0")
Maximum Load per Shelf		4 kg/shelf (8.8 lbs/shelf)	11 kg/shelf (24.3 lbs/shelf)	15 kg/shelf (33.1 lbs/shelf)
Electrical Configuration 110-130 VAC, 50/60 Hz	Nominal Power at 37 °C	10.6 W	42.9 W	49.5 W
	Max. Power Consumption	675.5 W	1184.3 W	1727.9 W
	Full Load Amps	5.3 A	9.2 A	13.4 A
Electrical Configuration 220-240 VAC, 50/60 Hz	Nominal Power at 37 °C	12.5 W	46.2 W	50.7 W
	Max. Power Consumption	598.8 W	1008.9 W	1270 W
	Full Load Amps	2.5 A	4.2 A	6.5 A
Contamination Control				
Contamination Control Methods		1) Main body is electrogalvanized steel with ISOCIDE™ antimicrobial coating; 2) HPA Validated 90°C Moist Heat Decontamination Cycle (Water in the external chamber must be drained first); 3) 0.2 micron inlet filter for gas inputs; 4) ULPA filter*****		

* Data recorded under optimum factory setting conditions.

** For temperature not exceeding 37°C.

*** For CO₂ not exceeding 5.2%.

**** Esco does not guarantee condensation-free chamber at humidity level higher than 90%.

***** Not available for 50 L.



CelMate®

CO₂ Incubators

INTRODUCTION

Esco now offers the new CelMate®, 170-liter and 240-liter, entry-level cell culture CO₂ incubator with superb contamination control.

This is specifically designed for laboratories looking for a cost-effective CO₂ incubator that can provide the best protection for their cell culture.



CelMate® CO₂ Incubators available in 170 L and 240 L.

CONTAMINATION CONTROL SYSTEMS:

- ULPA Filtration System
- Validated SwiftCon™ 90°C Moist Heat Decontamination Cycle (20 hours)
- Filtered Gas Injection Lines
- ISOCIDE™ Antimicrobial Exterior Coating

OPTIONS AND ACCESSORIES:

All options and accessories for CelCulture® CO₂ incubators are also appropriate to use on CelMate® CO₂ incubator.

ORDERING INFORMATION

MODELS	ITEM CODE	DESCRIPTION
CLM-170B-8	2170106	CelMate® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-8	2170107	CelMate® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CLM-170B-9	2170250	CelMate® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-9	2170251	CelMate® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)

GENERAL SPECIFICATIONS

CELMATE® CO₂ INCUBATORS

GENERAL SPECIFICATIONS		CLM-170B-__		CLM-240B-__	
TEMPERATURE					
Temperature Control Method		Direct Heat and Air Jacket using Microcontroller PI			
Ambient Temperature Range		18 to 34°C (64 to 93 °F)			
Temperature Range, °C		Ambient +3 to 60			
Temperature Uniformity*, °C		< ±0.5			
Temperature Accuracy, °C		<± 0.1			
Recovery Time** (after 1 min. door opening, 97% from initial value)		≤5 minutes			
CO ₂					
CO ₂ Control System		Microcontroller PI			
CO ₂ Range, % CO ₂		0-20			
CO ₂ Accuracy, % CO ₂		±0.1			
CO ₂ Sensor		Infrared (IR) Sensor			
CO ₂ Recovery Time*** (after 1 min. door opening)		≤4 minutes			
HUMIDITY					
Humidification Method		Humidity pan			
Humidity Range****, % RH		Up to 95%			
PHYSICAL CONSTRUCTION					
Interior Volume		170 L (6 ft³)		240 L (8.5 ft³)	
External Dimensions (W x D x H)		660 x 660 x 900 mm (26.0" x 26.0" x 35.4")		750 x 770 x 900 mm (29.5" x 30.3" x 35.4")	
Internal Dimensions (W x D x H)		505 x 535 x 633 mm (19.9" x 21.1" x 24.9")		595 x 640 x 633 mm (23.4" x 25.2" x 24.9")	
Chamber Construction	Main Body	Electrogalvanized steel with ISOCIDE™ antimicrobial coating			
	Interior Material	Stainless steel, type 304			
	Number of Shelves	4			
	Maximum Number of Shelves	7			
	Shelves Area (W x D)	465 x 470 mm (18.3" x 18.5")		550 x 560 mm (21.7" x 22.0")	
	Maximum Load per Shelf	11 kg/shelf (24.3 lbs/shelf)		15 kg/shelf (33.1 lbs/shelf)	
Electrical Configuration 110-130 VAC, 50/60 Hz	Nominal Power at 37°C	42.9 W		49.5 W	
	Maximum Power Consumption	1184.3 W		1727.9 W	
	Full Load Amps	9.2 A		13.4 A	
Electrical Configuration 220-240 VAC, 50/60 Hz	Nominal Power at 37°C	46.2 W		50.7 W	
	Maximum Power Consumption	1008.9 W		1270 W	
	Full Load Amps	4.2 A		6.5 A	
Shipping Weight		120 kg (264.6 lbs)		155 kg (341.7 lbs)	
Shipping Dimensions (W x D x H)		850 x 720 x 1120 mm (33.5" x 28.3" x 44.1")		850 x 850 x 1120 mm (33.5" x 33.5" x 44.1")	
Shipping Volume		0.70 m³ (24.85 ft³)		0.79 m³ (28.03 ft³)	
CONTAMINATION CONTROL					
Contamination Control Methods		1) Main body is electrogalvanized steel with ISOCIDE™ antimicrobial coating; 2) 90°C Moist Heat Decontamination Cycle (no decon pump); 3) 0.2 micron inlet filter for gas inputs; 4) ULPA filter			

*Data recorded under optimum factory setting conditions.

**For temperature not exceeding 37°C .

***For CO₂ not exceeding 5.2%

**** Esco does not guarantee condensation-free chamber at humidity level higher than 90%.

CelSafe®: NEW GENERATION CO₂ INCUBATOR

Esco's CelSafe® CO₂ incubator with touch screen user interface and latest advanced technology represents safety of your precious samples, efficiency on your lab work and enhanced user experience.

With CelSafe®, you will never look for another CO₂ / O₂ incubator




KEY COMPONENTS OF CelSafe® CO₂ INCUBATOR




MAGNETIC DOOR LOCK

- Manual and automatic lock functions.



STERILIZATION COOLING FAN

- Helps to emit the hot air during sterilization cycle.



DIRECT HEAT AND AIR JACKETED DESIGN

- Fast uniform heating
- Rapid temperature recovery without overshoot
- Air jacket improves chamber stability
- Double insulation system



INNER DOOR LATCH

- Locks the glass door
- Automatically turns off the pump, gas supply, and heating functions when inner door is opened



SHELVING

- Perforated shelving to improve uniformity
- Anti-tip
- Stainless steel
- Built-in grip



WATER RESERVOIR (Active Humidification Mode)

- User can precisely set %RH required for specialized application.




WATER PAN (Standard Model)

- Precisely heated by base heater to provide humidity



LEVELING FEET

- Easily adjustable






PRECISE PARAMETER CONTROL

HOMOGENEOUS CHAMBER ENVIRONMENT


- Direct heat temperature control system
- Excellent temperature uniformity under normal operating conditions
- Gentle forced convection design for even parameter distribution



LATEST INFRARED CO₂ SENSOR TECHNOLOGY

The new carbon dioxide IR Sensor probe withstands high temperature sterilization.

- CARBOCAP® technology for heat durability and long term stability
- With full temperature and pressure compensation
- No need to remove CO₂ probe inside the chamber during sterilization cycle, thus saving time and reducing risk of cross contamination.
- Water vapor, pressure, dust, other chemicals and gases, and temperature changes cannot affect the performance of the IR sensor.
- Internal pressure sensor improves accuracy and stability.
- Sensor head is heated to prevent condensation.



180°C HIGH HEAT STERILIZATION SYSTEM

With a simple touch on the screen, CelSafe® sterilization cycle assures deactivation of microbes, spores, fungi, vegetative cells and other harmful microorganisms that can affect the growth of your precious samples.

MICROORGANISM	Before Decon	After Decon
<i>Bacillus atrophaeus</i>	Log 6	0
<i>Aspergillus brasiliensis</i>	Log 4	0
<i>Pseudomonas aeruginosa</i>	Log 6	0
<i>Staphylococcus epidermidis</i>	Log 6	0
<i>Escherichia coli</i>	Log 6	0
<i>Staphylococcus aureus</i>	Log 6	0
<i>Enterococcus faecalis</i>	Log 6	0

- Fully automatic high heat sterilization cycle with a simple touch on the screen.
- Effectiveness of high heat sterilization cycle is validated through in-house laboratory test.
- All components and accessories are designed to meet high temperature requirement.
- Complies with different international guideline requirements for dry heat sterilization such as U.S. and E.U. Pharmacopeias.
- Everything is standard to avoid running cost on other external accessories and consumables to perform sterilization cycle.



OPTIMIZED CLEAN CHAMBER DESIGN

Less components mean more space for your samples.

- New ribbed design chamber allows installation of shelves without screws or pilasters.
- Minimize risk of contamination
- Easy maintenance
- Quick and easy to clean
- More chamber space

MODELS	CelCulture®	CelSafe®
90 mm Petri Dish	480 pcs	700 pcs
Treated Flask 25 cm ² Surface Area	384 pcs	384 pcs
Treated Flask 175 cm ² Surface Area	48 pcs	64 pcs
Cell Culture Plate (96 wells)	240 pcs	300 pcs
Cell Culture Plate (24 wells)	240 pcs	300 pcs
Cell Culture Plate (48 wells)	240 pcs	300 pcs




USER, SAMPLE, AND ENVIRONMENT PROTECTION

ISOCIDE™ ANTIMICROBIAL POWDER COATING


- Main body is electro-galvanized steel with white oven-baked epoxy-polyester antimicrobial powder coating.
- Isocide™ cannot be washed out or diminished by repeated cleaning.
- Ensures healthier, safer and cleaner lab environment.

ISOCIDE™ Antimicrobial Coating Performance




GREEN PRODUCT

DOUBLE INSULATION SYSTEM = LESS HEAT EMISSION




Regular CO₂ Incubator with Single Insulation
 Heat Emission Value: Approximately 42 W/sec
 Energy Consumption: 150 kWh




CelSafe® CO₂ Incubator with Double Insulation
 Heat Emission Value: Approximately 39 W/sec
 Energy Consumption: 142 kWh

COMPLETE SECURITY SYSTEM

- Multiple over-temperature protection system guarantees maximum sample, user and environment protection.
- All electrical components are UL-recognized.
- Electrical circuit protection is in accordance with UL requirements.
- PIN code prevents unauthorized access on screen menu and functions.
- Magnetic Door Lock System
 - Manually locks during normal operation to protect samples
 - Automatically locks during high heat sterilization cycle to protect users
- Inner Door Latch function turns off pump, gas supply, and heating functions when the door is opened.



Door Lock Option




PIN Code Security Display

CELTOUCH™ SCREEN CONTROL SYSTEM

High-tech, Simple and Functional CelTouch™ screen interface


- Big, clear, and easy-to-read parameter display
- Easy to follow on-screen icon menus
- Actual Data Graph, Data Logging functions, Event Logs and Alarm Functions are easily seen on the screen.
- Easy download of data log files using USB Write menu
- Fully functional even with gloved fingers
- Multiple language selections: English, German, Spanish, French, Italian



Home Screen



Icon Menu



Graph Display

The screen displays the 'ESCO' logo at the top. Below it, there is a table with the following data:

Date and Time	T °C	%CO ₂	%O ₂	%RH
20 April 2015 - 11:40	37.0	5.0	5.1	85.0
20 April 2015 - 11:39	37.0	5.0	5.1	85.0
20 April 2015 - 11:38	37.0	5.0	5.1	84.0
20 April 2015 - 11:37	36.9	5.1	5.0	84.8
20 April 2015 - 11:36	36.9	5.0	4.9	84.9
20 April 2015 - 11:35	37.0	5.0	5.0	85.0
20 April 2015 - 11:34	37.0	5.0	5.0	85.0
20 April 2015 - 11:33	37.0	5.0	5.0	85.0

At the bottom, there is a navigation bar with icons for home, menu, add, help, and speaker, along with the date and time: Wed, 22 April 2015 14:00.

Data Log Display



Language Option Display

COMPLETE DATA COLLECTION WITH USB WRITE FUNCTION

Easy Download of Data and Event Logs with CelTouch™ USB Write Feature

USB PORT



USB Port




USB Write Screen

ANALOG OUTPUT

- Stand-by 0-5 VDC 4-20 mA analog output allows the chamber to be connected to an in-house data acquisition or alarm system.

ALARM CONTACTS

- A set of relay contacts located on the rear of the unit is provided to monitor temperature, humidity or CO₂ alarms. The alarm contacts can be connected to a remote alarm system.



Analog Output and Alarm Contacts

FLEXIBILITY ON YOUR CELSAFE® CO₂ INCUBATOR

ACTIVE HUMIDIFICATION SYSTEM

- In order to provide optimal environmental conditions for cell growth that requires specific relative humidity, the CelSafe® CO₂ incubator with optional active humidity control allows user to actively control humidity from up to 95%.
- Water reservoir is located at the back of the chamber.
- Heats up and maintains water reservoir based on RH control.
- Water inlet valve is triggered by water level sensor.



Water Reservoir

SUPPRESSED O₂ MODEL

- Esco CelSafe® CO₂ incubators with suppressed O₂ provide accurate environmental control inside the incubator chamber. Oxygen levels are controlled through precise introduction of nitrogen into the incubator culturing system.
- Esco CelSafe® CO₂ incubator with suppressed O₂ has an optional nitrogen gas switching system in order to install two nitrogen gas tanks to make sure you will not run out of N₂ gas.
- New Zirconium type O₂ Sensor provides faster response time and greater reliability than Galvanic type O₂ Sensor.

STILL WANT MORE PROTECTION?

ESCO GOT YOU COVERED USING ESCO PROTECT® SOFTWARE SYSTEM


PROtect®

A completely independent and redundant sample monitoring system, which is a critical component in providing protection for important sample.

CFR-21 Compliant Software

WHAT ESCO PROTECT® DOES?

Esco PROtect® monitoring system includes radio-frequency recorders which can be connected to any brands of freezer, cold-storage system and temperature-controlled equipment. These recorders automatically send data to a central server to monitor temperature and other parameters.



LEGEND:



ESCO Protect
UN/DG



ESCO Protect
ALARM



Visual and
Audible Alert



ESCO Protect
EMD-1





Email Alert

Ethernet network

Radio Communication
(200-300 meters)

REAR PANEL



1 Cooling Fan

The cooling fan prevents the electrical panel from overheating.




2 Power Supply Inlet

The power supply inlet connects the incubator unit to the power source.



3 ModBus LAN Port

The ModBus LAN Port provides serial communication port for PC. It can be daisy-chained from product to product and connected to a PC.



4 Analog Port

The analog port allows the incubator to output analog signals representing temperature, CO₂/O₂ concentration and relative humidity, depending on the options available in the incubator. This allows the incubator to be connected to an in-house data acquisition or alarm system.



5 Alarm Contact

A set of relay contacts located on the rear of the unit is provided to monitor temperature, humidity or CO₂ alarms. The alarm contacts can be connected to a remote alarm system.



6 N₂ Gas Supply Inlet (for Suppressed O₂ model)

The N₂ gas supply inlet is only applicable for models with N₂ control function. Inlet pressure requirement is 15 psi.



7 CO₂ Gas Supply Inlet

The CO₂ gas supply inlet connects the CO₂ gas supply with the incubator unit. Inlet pressure requirement is 15 psi.



8 Gas Inlet Filter

In-line filters are provided to remove any contaminants from the gas supply.



9 Access Port

Allows cables, hoses or additional sensors to be routed into the work space. A rubber stopper with controlled leak is installed as standard configuration and is part of standard accessories.

ORDERING INFORMATION

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CLS-170B-8	2170187	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ Control, High Heat Sterilization, 230 VAC, 50/60 Hz
CLS-170B-9	2170188	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ Control, High Heat Sterilization, 115 VAC, 50/60 Hz

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER AND ACTIVE HUMIDIFICATION

MODELS	ITEM CODE	DESCRIPTION
CLS-170B-8-RH	2170192	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ Control, High Heat Sterilization, Active Humidification, 230 VAC, 50/60 Hz
CLS-170B-9-RH	2170194	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ Control, High Heat Sterilization, Active Humidification, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER

MODELS	ITEM CODE	DESCRIPTION
CLS-170T-8	2170130	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ and O ₂ Control, High Heat Sterilization, 230 VAC, 50/60 Hz
CLS-170T-9	2170151	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ and O ₂ Control, High Heat Sterilization, 115 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER AND ACTIVE HUMIDIFICATION

MODELS	ITEM CODE	DESCRIPTION
CLS-170T-8-RH	2170193	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ and O ₂ Control, High Heat Sterilization, Active Humidification, 230 VAC, 50/60 Hz
CLS-170T-9-RH	2170195	CelSafe® Incubator, 170 L, IR Sensor, CO ₂ and O ₂ Control, High Heat Sterilization, Active Humidification, 115 VAC, 50/60 Hz

GENERAL SPECIFICATIONS	
CelSafe® CO ₂ INCUBATORS	
CLS-170_ _	
TEMPERATURE	
Ambient Temperature Range	18°C to 34°C (64°F to 93 °F)
Temperature Control Method	Direct Heat & Air Jacket using Microcontroller PI
Temperature Control Range, °C	Ambient +3 to 60
Temperature Uniformity*, °C	±0.5
Temperature Accuracy, °C	0.1
CO ₂	
CO ₂ Control System	Microcontroller PI
CO ₂ Control Range, % CO ₂	0-20
CO ₂ Accuracy, % CO ₂	0.1
CO ₂ Fluctuation, % CO ₂	± 0.2
CO ₂ Sensor	IR sensor
O ₂ SPECS (FOR SUPPRESSED O ₂ MODEL)	
O ₂ Control System	Microcontroller PI
O ₂ Control Range, % O ₂	1-20.7%
O ₂ Accuracy, % O ₂	0.1
O ₂ Fluctuation, % O ₂	± 0.2
O ₂ Sensor	Zirconium Oxide
HUMIDITY	
Humidification Method	Humidity pan (Natural Humidification) / Humidity reservoir (Optional Active Humidification)
Humidity Range**	Standard Unit: 85% - 90% Active Humidification (Optional): 90% - 95%
PHYSICAL CONSTRUCTION	
Interior Volume	170 L (6.0 ft ³)
External Dimensions (W x D x H)	660 x 720 x 980 mm (26.0" x 28.3" x 38.6")
Internal Dimensions (W x D x H)	505 x 530 x 635 mm (19.9" x 20.9" x 25.0")
Net Weight	99 Kg (218 lbs)
Shipping Weight	118 Kg (260 lbs)
Shipping Dimensions (W x D x H)	850 x 770 x 1135 mm (33.5" x 30.3" x 44.7")
Exterior Material	Electrogalvanized steel with ISOCIDE™ antimicrobial coating
Interior Material	Stainless steel, type 304
Number of Shelves	4
Maximum No. of Shelves	7
Shelves Area (W x D)	502 mm x 475 mm (19.8" x 18.7")
Max. Load per Shelf	11 Kg/shelf (24.3 lbs/shelf)
Available Electrical Configuration	230 VAC, 50/60 Hz (8)
	115 VAC, 50/60 Hz (9)
Interior Material	Stainless Steel, Type 304
CONTAMINATION CONTROL	
Contamination Control Methods	1) 180°C High Heat Sterilization Cycle; 2) Main body is electrogalvanized steel with ISOCIDE™ antimicrobial coating; 3) 0.2 micron inlet filter for gas input; 4) 1 micron in-line filter for air circulation

* Data recorded under optimum factory setting conditions

** Esco does not guarantee condensation-free chamber at humidity level higher than 90%.

OPTIONS AND ACCESSORIES




COA-1001 / COA-1001-F Humidity Display

This option allows the incubator to monitor the relative humidity inside the chamber. The probe for the sensor works in freezing conditions (-70°C) and also in temperatures up to 180°C . The sensor is easy to install and has excellent accuracy. The airflow in the chamber does not affect the measurement. The sensor is maintenance-free. It does not need to be removed during 90°C moist heat decontamination cycle.



COA-1002 / COA-1002-F CO₂ Backup

This option allows two tanks of CO₂ to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low gas pressure is detected on the primary tank.



COA-1005 / COA-1005-F Analog Output

A set of relay contacts is provided at the rear of the incubator that allows the incubator to output analog signals representing the temperature, CO₂ / O₂ content and relative humidity, depending on the options available in your incubator. This allows the chamber to be connected to an in-house data acquisition or alarm system. This option can also be field-installed.

The analog signal outputs can be set to operate in either voltage DC (0-5 Vdc) or current (4-20 mA) mode. The factory default setting is voltage. Switch on the board to toggle between the modes.



COA-2030/ COA-2030-F Sealed Inner Door Kit with 2 glass doors (50L) COA-1006/ COA-1006-F Sealed Inner Door Kit with 4 glass doors (170L) COA-2029/ COA-2029-F Sealed Inner Door Kit with 4 glass doors (240L) COA-2040/ COA-2040-F Sealed Inner Door Kit with 6 glass doors (240L) COA-2038 / COA-2033-F Sealed Inner Door Kit for CelSafe®

CelCulture® CO₂ incubators can be equipped with 2, 4 or 6 glass doors, that can be opened horizontally which allows access to defined sections of the incubator without affecting much the inner atmosphere of the chamber. This minimizes recovery time and contamination risks. The sealed-inner door is also reversible as same as the outer door which can be installed to be opened either from-right-to-left or from-left-to-right. The sealed-inner door is available as a factory-installed option or field installed retrofit kit.



COA-1007 / COA-1007-F N₂ Back-up

This option allows two tanks of N₂ to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low gas pressure is detected on the primary tank.



COA-2018-F (50L) / COA-2001-F (170 L) / COA-2019-F (240 L) / COA-2034-F (CLS) Roller Base

Roller base is available with casters for mobility of your incubators and to provide protection against floor contamination.



COA-2020-F (50L) / COA-2002-F (170 L) / COA-2021-F (240 L) / COA-2035-F (CLS) Floor Stand 200 mm (8.0") With Adjustable Feet

Floor stands are available with adjustable feet, with a nominal range of 180 mm to 250 mm (7.1" to 9.8") for comfortable access to the incubator and to avoid floor contamination.



COA-2022-F (50L) / COA-2003-F (170 L) / COA-2023-F (240 L) / COA-2036-F (CLS) Floor Stand 700 mm (27.6") with Casters

This support stand raises the incubator to a height of 700 mm (27.6") above the floor for comfortable access. It comes with casters for mobility of your incubators.



COA-2005-F 2-Stage Gas Regulator for CO₂/N₂

CO₂ and N₂ gas input regulators reduce pressure from the tank to the incubator. It has dual pressure gauges, barbed line connection and shut-off valve. It prevents over-pressurization of the gas supply into the incubator which could cause the tubing to burst.

- CGA 320 connector (U.S. Standard)
- BP-BS341-#8-NT4 connector (British Standard)
- G5/8-RH connector (China Standard)

Note: Compatible with European DIN477, French NFE29-650 and Australia AS2473



COA-2024-F (50L)/ COA-2007-F (170 L)/ COA-2025-F (240 L)/ COA-2037-F (CLS) Extra Shelf (Stainless Steel) for Standard Stainless Steel Chamber

COA-2026-F (50L) / COA-2027-F (170 L) / COA-2028-F (240 L) Extra Shelf (Copper) for Standard Copper Chamber

Each CelCulture® CO₂ incubator comes standard with 3 shelves for 50 L / 4 shelves for 170 L & 240L and it can accommodate up to a maximum of 4 shelves for 50 L / 7 shelves for 170 L & 240 L.



COA-2008-F Stacking Kit

The stacking kit is a provision to stack one incubator on top of another incubator. Four stacking brackets are included as standard inside the Accessories Kit Box with each incubator.



COA-2010-F Electronic CO₂ Analyzer, For CO₂ / Temp Measurement

COA-2016-F Electronic CO₂ + O₂ Analyzer, For CO₂ / O₂ / Temp Measurement

COA-2017-F Electronic CO₂ + O₂ + RH Analyzer, For CO₂ / O₂ / RH / Temp Measurement

The electronic analyzer allows the measurement of CO₂ concentration, O₂ concentration, relative humidity and temperature (temperature probe already included).




COA-2012-F 6" Chart Recorder, Temp, 115/230 VAC, 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature data.



COA-2013-F 8" Chart Recorder, Temp/Temp, 115/230 VAC, 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 8" chart of temperature data and comes with 2 remote probes for dual temperature monitoring.



COA-2014-F 6" Chart Recorder, Temp/RH, 115/230VAC 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature and humidity data.



COA-2015-F Inner Door Shelving Kit


(4 Sets With Total 12 Mini-Shelves For One Incubator) (170 L)

These mini-shelves are to be used with the Sealed Inner Door Kit installed. There are 4 sets with a total of 12 mini-shelves on each incubator.



5250001 Voyager® Software Kit

Esco Voyager® is a PC-based software package developed for the remote monitoring, datalogging and programming / device configuration of Esco controlled environment laboratory equipment. Compatible equipment includes laboratory ovens and incubators, low temperature incubators, CO₂ incubators, and ultra-low temperature freezers.



COA-2004 CCL-170 / COA-2039 CLS-170 2-UNITS FLOOR STAND

This floor stand allows two units to be stacked without being physically in contact with each other. For the lower unit, it uses roller base for mobility and for easy pull out of the lower unit without the need to remove the upper unit in case of troubleshooting. Floor stand for upper unit also has casters for easy relocation.

ORDERING INFORMATION

ACCESSORIES	ITEM CODE	DESCRIPTION
COA-1001	5170470	Humidity Display, Factory-installed
COA-1001-F	5170471	Humidity Display, Field-installed Kit
COA-1002	5170472	CO ₂ Backup (Tank Switcher), Factory-installed
COA-1002-F	5170473	CO ₂ Backup (Tank Switcher), Field-installed
COA-1004	5170474	Reversed Door Swing, Factory-installed
COA-1005	5170475	Analog Outputs, Factory-installed
COA-1005-F	5170476	Analog Outputs, Field-installed
COA-2030	5170672	Sealed Inner Door Kit for 50 L (2 Glass Doors with Latches), Factory-installed
COA-2030-F	5170673	Sealed Inner Door Kit for 50 L (2 Glass Doors with Latches), Field-installed
COA-1006	5170477	Sealed Inner Door Kit for 170 L (4 Glass Doors with Latches), Factory-installed
COA-1006-F	5170488	Sealed Inner Door Kit for 170 L (4 Glass Doors with Latches), Field-installed
COA-2029	5170654	Sealed Inner Door Kit for 240 L (4 Glass Doors with Latches), Factory-installed
COA-2029-F	5170655	Sealed Inner Door Kit for 240 L (4 Glass Doors with Latches), Field-installed
COA-2040	5170783	Sealed Inner Door Kit for 240 L (6 Glass Doors with Latches), Factory-installed
COA-2040-F	5170786	Sealed Inner Door Kit for 240 L (6 Glass Doors with Latches), Field-installed
COA-2038	5170701	Sealed Inner Door Kit for CLS (4 Glass Doors With Latches), Factory Installed
COA-2033-F	5170696	Sealed Inner Door Kit for CLS (4 Glass Doors With Latches), Field Installed
COA-1007	5170490	N ₂ Back-up (Tank Switcher), Factory-installed
COA-1007-F	5170491	N ₂ Back-up (Tank Switcher), Field-installed
COA-2018-F	5170419	Roller Base (50 L)
COA-2001-F	5170478	Roller Base (170 L)
COA-2019-F	5170420	Roller Base (240 L)
COA-2034-F	5170697	Roller Base (CLS)
COA-2020-F	5170421	Floor Stand 200 mm (8.0") with Adjustable Feet (50 L)
COA-2002-F	5170479	Floor Stand 200 mm (8.0") with Adjustable Feet (170 L)

COA-2021-F	5170422	Floor Stand 200 mm (8.0") with Adjustable Feet (240 L)
COA-2035-F	5170698	Floor Stand 200 mm (8") Adjustable Feet (CLS)
COA-2022-F	5170423	Floor Stand 700 mm (27.6") with Casters (50 L)
COA-2003-F	5170480	Floor Stand 700 mm (27.6") with Casters (170 L)
COA-2023-F	5170424	Floor Stand 700 mm (27.6") with Casters (240 L)
COA-2036-F	5170699	Floor Stand 700 mm (27.6") with Casters (CLS)
COA-2005-F	5170481	2-Stage Gas Regulator for CO ₂ / N ₂ Choose one of the connectors below: 1080588 - CGA 320 Connector (US standard) 1080589 - BP-BS34-#8-NT4 Connector (British standard) 1080590 - G5/8-RH Connector (China standard)
COA-2024-F	5170425	Extra Shelf (50 L, Stainless Steel)
COA-2007-F	5170327	Extra Shelf (170 L, Stainless Steel)
COA-2025-F	5170426	Extra Shelf (240 L, Stainless Steel)
COA-2037-F	5170700	Extra Shelf (CLS, Stainless Steel)
COA-2026	5170427	Extra Shelf (50 L, Copper)
COA-2027	5170428	Extra Shelf (170 L, Copper)
COA-2028	5170495	Extra Shelf (240 L, Copper)
COA-2008-F	5170483	Stacking Kit (one set included with every unit purchased)
COA-2010-F	5170329	Electronic CO ₂ Analyzer, For CO ₂ / Temp Measurement (with Temperature Probe)
COA-2016-F	5170397	Electronic CO ₂ + O ₂ Analyzer, For CO ₂ / O ₂ / Temperature Measurement
COA-2017-F	5170398	Electronic CO ₂ + O ₂ + RH Analyzer, For CO ₂ / O ₂ / RH / Temperature Measurement
COA-2011-F	2170020	IQ / OQ Documentation
COA-2012-F	2170021	6" Chart Recorder, Temp, 115/230 VAC, 50/60 Hz
COA-2013-F	2170022	8" Chart Recorder, Temp, 115/230 VAC, 50/60 Hz
COA-2014-F	2170023	6" Chart Recorder, Temp/RH, 115/230 VAC, 50/60 Hz
COA-2015-F	5170487	Inner Door Shelving Kit for 170 L (4 sets with total 12 mini-shelves for one incubator)
COA-2004-F	5170489	2-units Floor Stand Stacking Kit (170L)
COA-2039-F	5170778	2-units Floor Stand Stacking Kit (CLS)
Voyager®	5250001	Voyager® Software Kit


Advanced Cell Culture with Esco CO₂ Incubator and CelCradle™

CELCRADLE™ : CRADLE FOR HIGH DENSITY CELLS

CelCradle™ is a cost-effective, single-use benchtop bioreactor system capable of supporting high density culture of adherent cells. It is designed based on the concept of bellow-induced intermittent flow of media and air through porous matrices, where cells reside. This provides a low shear stress, high aeration, and foam-free culture environment.

During operation, the CelCradle™ bottle is partially filled with media and inoculated with cells. The media is raised and lowered alternately to submerge and expose the matrices, creating a dynamic interface between air and media on cell surface to maximize nutrient uptake and oxygen transfer.

CelCradle™ system is part of the tide motion bioreactor system, which features linear scalability up to 5,000 L packed bed volume!



BioNOC™ II: Heart of the Tide Motion System


BioNOC™ II are macrocarriers that allow attachment of cells for adherent cell culture. Each CelCradle™ bottle consists of 5.5g of BioNOC™ II carriers, providing 15,000 cm² of surface area for attachment and growth of cells. Apart from its high surface area, BioNOC™ II carriers feature enhanced biocompatibility, long hydrophilicity, high porosity, low lint content, and excellent mechanical strength. The characteristics of BioNOC™ II, combined with the tide motion principle, allow the CelCradle™ system to support the high density culture of adherent cells.

Features:

- Stainless Steel 304 L BA CelCradle™ Stage capable of operating 4 CelCradle™ bottles simultaneously and compatible with a CO₂ incubator
- Pre-sterilized and ready-to-use disposable CelCradle™ bottles
- Provides a low shear stress and foam-free culture environment that has no O₂ limitation
- A single CelCradle™ bottle has the same productivity of up to 20 850 cm² roller bottles
- Compact design allows the CelCradle™ to be placed inside a 6 ft³ CO₂ incubator
- Easy parameter optimization
- Capable of performing batch, fed-batch or perfusion culture mode
- Compatible with most media formulations including serum-free media
- Contains BioNOC™ II carriers with specially treated surface to allow growth of most anchorage-dependent cells. BioNOC™ II also allows easy harvest of whole cells, cell components or secreted proteins
- Easy-scale up by using additional bottles or by using TideCell® bioreactor system

Applications:

- Human and Animal Vaccines
- Autologous and Allogeneic Cell Therapy
- Culture of anchorage-dependent/ adherent cells
- Overcome limitations of stirred tank bioreactors using microcarriers
- Conversion from Roller Bottles to closed system, single-use cell culture
- Mammalian and insect cell research
- Monoclonal antibody production
- Protein production



CELCRADLE™ BATCH SYSTEM

Useful for batch and semi-batch operation where process components are easily traceable.

Simple to operate - virtually no learning curve

Controller - adjusts nutrient and gas exchange of the bottle

Magnetized controller enables convenient positioning on the outside surface of the incubator



Bottle cap is integrated with 0.22µm ventilation filter

BioNOC™ II - acts as the packed bed where cells can adhere to, simplifying media replacement and product harvesting

Sampling Port - enables aseptic removal of BioNOC™ II macrocarriers for cell counting

Retaining Ring - locks bottle in place in the CelCradle Stage

Collapsible bellows

CELCRADLE™ CONTINUOUS SYSTEM

Extremely compact system fits in most standard CO₂ incubators.

Easy to read LED display

Up: 1.5 mm/s
T_H: 04 M 04 S ^



Useful for continuous operation where process components are easily traceable

Tubing set with peristaltic pump head enable recirculation and continuous feeding of nutrient

CelFeeder enables 4-pump operation with individual programming setting

Autoclavable pump head supports consistent, calibration-free feeding rate


DISPOSABLE CELCRADLE™ BOTTLES

A complete product line of CelCradle™ meets your specific needs. Different CelCradle™ bottles cover 90% of applications in cell culture.

- Batch, semi-batch or continuous culture
- BioNOC™ II carriers or preferred microcarriers
- Cell harvest with or without trypsin
- Disposable Bottle

Bottle	Item Code	Secreted Protein, viruses (adherent cells)	Cell Harvest (for non-secreted proteins, viruses or cells)	Carrier Harvest (for protein extraction or reuse of carriers)
CelCradle™ 500	1400001	Best Application	Applicable, but not optimal	Applicable, but not optimal
CelCradle™ 500A	1400003	Applicable, but not optimal	Best Application	Best Application
CelCradle™ 500P	1400002	Best Application	Applicable, but not optimal	Applicable, but not optimal
CelCradle™ 500AP	1400004	Applicable, but not optimal	Best Application	Best Application

OPTIONS AND ACCESSORIES



CelFeeder

The CelFeeder pump module is an auxiliary peristaltic pump used for the recirculation or perfusion processes for CelCradle™ 500 high density continuous cell culture system.



Tubing Complete Set

The Tubing Complete Set includes pre-assembled tubes, reusable pump head and head plate with a sampling port to support the continuous culture in CelCradle™-500P system.



Disposable Tubing Accessory

The Disposable Tubing Accessory provides simple options to replace the tubes in the Tubing Complete Set. This replacement avoids tubing deterioration during operation. It is recommended to replace the tubes after three (3) times of use.



GlucCell® Glucose Monitoring System

The GlucCell® enables simple and accurate glucose measurements using disposable test strips.



Crystal Violet Dye Nucleus Count Kit

The Crystal Violet Dye Nucleus Count Kit contains crystal violet dye, citric acid and detergent used to disrupt the cells and release cell nuclei for cell count. The CVD kit is an efficient reagent for cell count in a porous matrix.

ORDERING INFORMATION

Product Name	Item Code	Package
CelCradle™ System Complete	2230006	1 x CelCradle™ Stage 1 x GlucCell® Glucose Monitoring System
CelCradle™ Continuous System Complete	2230007	1 x CelCradle™ Stage 1 x GlucCell® Glucose Monitoring System 1 x CelFeeder Pump 2 x Tubing Complete Set
CelCradle™ Stage	2230005	1 x Main Console 1 x Control Box 1 x 100-240 V power adapter 1 x Signal Cable 1 x Manual CD 2 x Forceps 1 x Crystal Violet Dye Nucleus Count Kit
CelFeeder Pump	1400067	1 x CelFeeder Pump
Tubing Complete Set	1400011	1 x Disposable Tubing Accessory 1 x Pump Head 1 x Head Plate
Disposable Tubing Accessory	1400013	5 x Disposable Tubing Accessory
Disposable Tubing Set & Pump Head	1400012	1 x Tubing Set 1 x Pump Head
GlucCell® Glucose Monitoring System	1400009	1 x GlucCell® Glucose meter 2 x Glucose Test Strip Bottles (2 x 25 pcs)
GlucCell® Glucose Test Strip	1400010	2 x Glucose Test Strip Bottles (2 x 25 pcs)
Crystal Violet Dye Nucleus Count Kit	1400014	1 x CVD Bottle (100ml/bt)
Filtered Cap	1400015	Cap for CelCradle™ Bottle (Pack of 6)
Non-Vented Cap	1400016	Cap for CelCradle™ AP/P Bottle (Pack of 8)
Forceps	1400017	Used for sampling of BioNOC™ II carriers

After Sales Services

Parts Availability

Whenever service is needed and parts are required, minimizing downtime is a critical objective. Statistical usage analysis helps Esco to predict parts life, permitting Esco to manage logistics and stage proper inventories around the world. The combination of predictive maintenance, historical data and geospecific proximity assures our customers that parts and labor are available whenever service is scheduled through the local sales organization.

Registration, Documentation and Instruction

Quality control at Esco extends from research and development through engineering, manufacturing, shipment, delivery and customer feedback. Esco maintains an aggressive program to encourage warranty card registration by mail, email or online submission so that we know where Esco products are located and how they are being used. Rest assured that all information disclosed from warranty registrations will be kept confidential. All Esco products include unique serial numbers for identification. Documentation for all performance tests is archived and maintained for customer reference.

Online Technical Information

Site preparation instructions are useful before product arrival and installation. Installation and start-up manuals, operation manuals and quick reference guides are available anytime from the Esco resources online. An interactive online LiveSupport™ concierge center accessible through the Esco website offers extended hours of operation. LiveSupport™ permits users to dialogue directly with Esco personnel.

NSF International Accreditations

The National Sanitation Foundation (NSF) International is an independent, non-profit organization that provides standards development, product certification, auditing, education and risk management for public health and the environment.

In line with Esco's commitment in providing world class services worldwide, Esco has a large contingent of NSF accredited certifiers which makes Esco not only an **Excellent Standards Company** but also an **Excellent Service Company**, which exemplifies Esco's collective quest of being an **Eternally Successful Company**.

The NSF mark is your assurance that the product complies with all the standard requirements, tested by one of the most respected independent certification organizations in existence today. NSF conducts periodic unannounced inspections and product testing to verify that the product continues to comply with the standard. It is valued by consumers, manufacturers, retailers and regulatory agencies worldwide.


References and Links

For more information, you can visit Esco at www.escoglobal.com

Whatever your needs are,
**we serve it to you
on a silver platter.**



ESCO GLOBAL NETWORK



- Global Offices
- Joint Ventures
- Licensee
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers



ART Equipment
 Biological Safety Cabinets
 CO₂ Incubators
 Compounding Pharmacy Equipment
 Containment / Pharma Products
 Ductless Fume Hoods
 Freeze Dryers
 Lab Animal Research Products
 Laboratory Centrifuges
 Laboratory Fume Hoods
 Laboratory Ovens and Incubators
 Laboratory Shakers
 Laminar Flow Clean Benches
 PCR Cabinets
 PCR Thermal Cyclers
 Powder Weighing Balance Enclosures
 Ultra-low Temperature Freezers

The Esco Group of Companies is a global life sciences tools provider with sales in over 100 countries. The group is active in lab equipment, pharma equipment and medical devices. Manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, service and marketing subsidiaries are located in 12 major markets including the US, UK, Singapore, Japan, China and India. Regional distribution centers are located in the US, UK, and Singapore.

Life Science • Chemical Research • Assisted Reproductive Technology (ART) • Pharmaceutical Equipment • General Equipment

ESCO

WORLD CLASS. WORLDWIDE.

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777
 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com
 www.escoglobal.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
 Tel: +1 215-441-9661 • Fax 484-698-7757
 eti.admin@escoglobal.com • www.escolifesciences.us

Esco Global Offices: Bangladesh | China | Denmark | Hong Kong | India | Indonesia | Italy | Japan | Lithuania
 Malaysia | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam